

SYNCHRONOUS VIDEO INSTRUCTIONS

30-SECOND SIT TO STAND TEST

Instructions

1. <Clinician> will tell you when to start & stop.
2. Set chair against the wall.
3. Sit in the chair with your feet on the floor and facing the camera.
4. Cross your arms on your chest.
5. Stand up from the chair all the way then sit all the way back down as fast as you can.
6. Repeat for 30 seconds.
7. <Clinician says> Ready, Start. ... Stop.

Notes

TIPS

- Stop if you have any discomfort or feel unsteady.

EQUIPMENT:

- Chair: Straight back, no wheels, no arms (e.g., kitchen chair)
- Nearby Wall

HELPER: Optional

SETTING: Space so the video includes the entire body

VIDEO INFO

- Use a tripod or prop up phone as needed.
- Make sure the entire body can be seen.

TOTAL TIME: <5 minutes

PURPOSE: To test leg strength and endurance.

YOUTUBE VIDEO: [Watch a YouTube video of the assessment.](#)

Developed by Hoenig H, Jesus T, Hanson C, Hall K, Myers K, Romero S, Alexander N, Allen K, Moo L.
Project: Measuring Physical Function With Video Telehealth Technology, 2024-2025. Funded by the
Office of Connected Care.